

THE BENEFITS OF CITRASATE[™]

CITRATE-BASED DIALYSATE IN DIALYSIS TREATMENT





Is dialysate going to be the next focus of dialysis treatment?

For patients with end stage renal disease, the aim of dialysis treatment is removing the uremic toxins on the one hand and **correcting the levels of the blood electrolytes** on the other. ¹



The **dialysis fluid** plays a key role since it is composed of the elements (bicarbonate and acid concentrate) that **normalize the electrolyte constitution** of the blood.² Despite such an important role, the effect of the dialysis fluid on the patient is often underestimated.

Acetate-based dialysates are commonly used in dialysis treatment.³ For patients susceptible to the side effect of acetate such as **hemodynamic instability, metabolic acidemia and arterial hypotension**,⁴ acetate was reduced and replaced with a more biocompatible compound, citrate.⁵

The **benefits of citrate-based dialysis** is suggested specially for patients with **hemodynamic instability** (cardiovascular diseases such as arrythmias), with **higher inflammation, and fragile patients** with nutritional problems.^{6,7}

Due to the anti-coagulation properties of citrate, citrate-based dialysate is a suitable support for patients with **coagulation issues**.⁸

Citrate-based dialysate with high levels of magnesium is beneficial for patients with hypoalbuminemia and nutritional issues.^{7,9} Dialysis patients are complex and have various clinical needs. A suitable dialysis fluid tailored to the specific needs of the patient is an important aspect of a quality dialysis treatment.

What is Citrasate[™] and what are the benefits of the citrate-based dialysate?

Citrasate contains biocompatible citric acid (0.8 mmol/L), and 0.3 mmol/L sodium acetate as an acidifying agent.

This allows for reduction in the concentration of acetic acid from 3 mmol/L in traditional formulation to 0.3mmol/L

0.3mmol/l

Benefits of citrate-based dialysate

Higher survival

33%

Citrate-based dialysate is safe for the patients and increases the long-term survival of patients to up to 33%.⁷

100 0.06 0.06 0.06 0.06 0.07 0.00 0

Fig 1: Mortality based on dialysis fluid. Patients treated with citrate-based dialysate have significantly higher survival than patients treated with acetate-based dialysate during a 25-month follow-up (p= 0.008). Retrospective, single-center study, N= 137.

Correction of acidosis

Metabolic acidosis is a common complication in patients with chronic kidney disease.¹⁰

BETT

SIIRVIVAI

As shown in Fig 2, **Citrasate**, has a **significantly lower acidemia** post-dialysis vs. treatment with acetate.¹¹

Excess of bicarbonate in dialysis can lead to **post-HD alkalemia** (increase in pH) which is associated with **increased mortality**.¹²

The use of citrate **decreases post HD alkalemia vs. acetate**.²



Fig 2. Treatment with Citrasate in 13 patients during 3 dialysis treatments. T0: baseline, T2: 2H and T4: at the end of treatment.

Compared to acetate-based dialysate, the use of Citrasate normalizes the acidemia in patients and decreases the post-HD alkalemia.^{2,11}

Higher biocompatibility and reduction of inflammation

Inflammation and oxidative stress are risk factors for cardiovascular-related mortality in dialysis patients.¹³

Treatment with citrate-based dialysate or Citrasate improves the inflammation by:

- Reducing IL-6 (Fig 3)¹⁴
- Reducing serum CRP levels (Fig 4)⁸
- Decreasing eosinophils (higher biocompatibility)²





Fig 3. The IL-6 levels at baseline, T1 (1.5 months) and T2 (3 months). The pair-wise differences between T2 and baseline is significantly higher in acetate vs. citrate containing dialysate. p= 0.029. Multicenter, randomized, controlled study, N=47.

Fig 4. Mean serum CRP levels is significantly reduced in patients treated with Citrasate. Cross-over single center study, N=22.

Citrate-based dialysate leads to 3-fold lower CRP levels benefiting patients with cardiovascular diseases.^{8,13-15}

Improving the nutritional state

Dialysis by citrate improves nutrition by:

- Increasing the **appetite**⁶
- increasing the **lean mass index** (LTI)²
- Increasing serum albumin levels in patients with hypoalbuminemia²

HIGHER LTI

Treatment with **Citrasate significantly increases the serum albumin levels by 2.3%** even after a short period of 8 months compared with acetate (Fig 5).⁸



Fig 5. Mean serum albumin levels is significantly increased in patients treated with Citrasate. Cross-over single center study, N=22.

Reducing intra-dialytic hypotension

intra-dialytic hypotension occurs in 25% of patients in dialysis¹⁶ and is a significant predictor of mortality.¹⁷

After more than 4.000 treatments, the use of citrate-based dialysate shows a clear significant **reduction of 23%** in intra-dialytic hypotension.²

Cost-effectiveness for the dialysis center

Assuming an average time of 15 min for managing a patient with hypotension, this reduction of 23% **can save 0.05-0.35 euro per treatment**.

**



Fig 6. Percentage of arterial hypotension in 4.416 sessions of hemodialysis. Dialysis with acetate-based fluid has a significantly higher % of hypotension compared with citrate-based dialysate. Prospective, multicenter, crossover, randomized study, N=56.

A suitable alternative to heparin

Heparin is widely used as an anti-coagulation strategy for dialysis patients.¹⁸

However, heparin is associated with side effects and is contraindicated for certain patients.¹⁹⁻²¹

For these patients, the reduced dose of heparin can be compensated by citrate thanks to its **anti-coagulation** properties.²²

How does citrate inhibit the coagulation cascade?

Citrate ions, due to their negative electric charge, **chelate the positively-charged calcium**. Calcium is an essential component for the coagulation cascade and necessary for the formation of thrombin, an enzyme that catalyzes the formation of fibrin in the clot formation process. Forming a complex with calcium will render it inaccessible for the clotting cascade.²¹

Clinical evidence

The study by Meijers et al., shows the successful **termination of 94%** of sessions in both hemodialysis and pre-dilution hemodiafiltration with a citrate-based dialysate (in combination with a low-thrombogenic dialyzer, Solacea[™]).²³

Other studies have shown reaching **comparable dialysis efficiency with citrate** while gradually lowering the heparin doses.²²

Kossmann et al., switched 31 patients from normal dialysate containing acetate (NCD) to Citrasate and gradually reduced the heparin dosage. Even with the reduction of heparin by 55%, all the treatments with Citrasate successfully completed without clotting problems and the episodes of prolonged bleeding were decreased.²⁴



Fig 7a, 7b. Reduction of heparin by 55% decreases the bleeding episodes of >15 minutes. N=31.

Reduction of heparin cost

By using <u>Citrasate</u>, the dose of heparin can be reduced by up to 50%²⁵ which can result in an estimated **cost** reduction of 0.46 and 0.72 euro per treatment for unfractionated heparin (UF) and low molecular weight heparin (LMWH), respectively. These prices were calculated based on the UF price of 1.15 euro for the average of 8.000 IU, and the LMWH price of 1.80 euro for the average of 300 IU.

Heparin type	Estimated average reduction cost per treatment for 30%-50% heparin reduction
UFH	0.46 euro
LMWH	0.72 euro



Higher dialysis efficiency

The blood-dialyzer interaction causes the **clotting** of the dialyzer and the reduction of dialysis efficiency. By using citrate-containing dialysis fluid, calcium will be chelated and thus less available for the coagulation cascade.²¹

Clinical studies have shown **increased clearance** of small and middle molecules with **citrate-based** dialysate owning to its superior anti-clotting properties compared to acetate.^{26,27}

Improving anemia

Anemia is a common comorbidity in dialysis patients and associated with **cardiovascular mortality** and decreased quality of life.²⁸

Treatment with citrate-based dialysate is shown to **increase the hemoglobin levels** in anemic patients with a **1.000-unit reduction in erythropoietin stimulating agent** (ESA).⁹

Lower cost of ESA

Based on the average reduction of 1.000 units in ESA when using a citrate-based dialysate, **a range of 0.43-1.70 euro cost reduction** can be estimated depending on the type of ESA used. This reduction is based on the price of ESA ranging from 1.3-5.4 euro for the average dose of 1.000 units of ESA per treatment.

0.43-1.70 € SAVING PER TREATMENT

Benefits of Citrasate with higher magnesium (Mg)

In dialysis patients with serum Mg levels above 2.1 mg/dL, a **20% increase in 2-year survival** has been observed (Fig 8).⁷

The increased mortality with **decreased Mg** can be related to the following **consequences** of hypomagnesemia in dialysis patients:



Fig 8. Patients with serum Mg levels > 2.1 mg/dL have a significantly higher survival than patients with <2.1 mg/dL during a 25-month follow-up (p=0.008). Retrospective, single-center study, N= 137.

In patients receiving dialysis, on the one hand, **Mg is cleared** through the dialysis treatment and on the other hand, the **compensation** for Mg through dietary intake is usually **limited by the specific recommended diets**.⁷

In vitro, the higher Mg doses of >1 mmol/L improves the oxidative stress in immune cells of dialysis patients.³⁷

The current patterns indicate that higher levels of Mg in the dialysate could potentially improve patients' outcomes.

Using citrate-based dialysate has several benefits

For patients

Dialysis patients are complex and have different medical needs. **Individualization of the dialysis fluid** to specific medical needs of the patients has to be considered as an ideal solution.³⁸ Dialysis with citrate-based dialysate offers several benefits for patients.

For healthcare staff

- 33% increased survival of patients⁷
- **23% reduction of hypotensive episodes**, leading to **6 hours saved** for a dialysis unit with 100 patients

For dialysis centers

- Up to **0.35 euro** cost-saving on reduction of hypotensive episodes
- Up to **0.72 euro** cost-saving on heparin
- Up to 1.70 euro cost-saving on reduction of ESA



Total financial advantage of 2.77 euro per treatment



Why Citrasate?

Increased clearances Improvement of inflammation Improvement of nutrition Less intra-dialytic hypotension Acid-base control with less pre-dialysis acidosis

References

- 1. Coresh J. Update on the Burden of CKD. J Am Soc Nephrol. 2017
- 2. de Sequera P, et al. Advantages of the use of citrate over acetate as a stabilizer in hemodialysis fluid: A randomized ABC-treat study. Nefrologia (Engl Ed). 2021
- 3. Kohn OF, et al. Dual-concentrate bicarbonate-based hemodialysis: know your buffers. Artif Organs. 2012
- 4. Dhondup T, et al. Electrolyte and Acid-Base Disorders in Chronic Kidney Disease and End-Stage Kidney Failure. Blood Purif. 2017
- 5. Pérez-García R, et al. Citrate dialysate does not induce oxidative stress or inflammation in vitro as compared to acetate dialysate. Nefrologia. 2017
- 6. Matsuyama K, et al. Acetate-free blood purification can impact improved nutritional status in hemodialysis patients. J Artif Organs. 2011
- 7. Pérez-García R, et al. Hypomagnesaemia in haemodialysis is associated with increased mortality risk: its relationship with dialysis fluid. Nefrologia (Engl Ed). 2020
- 8. Weiss w, et al. SP681, Decrease of inflammation and bleeding times in chronic dialysis patients by long-term use of citrate-enriched dialysate. nephrology dialysis transplantation. 2017
- 9. Kuragano T, et al. Effects of acetate-free citrate-containing dialysate on metabolic acidosis, anemia, and malnutrition in hemodialysis patients. Artif Organs. 2012
- 10. Rezende LR, et al. Metabolic acidosis in hemodialysis patients: a review. J Bras Nefrol. 2017
- 11. Panichi V, et al. Post-dilution on line haemodiafiltration with citrate dialysate: first clinical experience in chronic dialysis patients. Scientific World Journal. 2013
- 12. Tentori F, et al. Association of dialysate bicarbonate concentration with mortality in the Dialysis Outcomes and Practice Patterns Study (DOPPS). Am J Kidney Dis. 2013
- 13. de Mutsert R, et al. Excess mortality due to interaction between protein-energy wasting, inflammation and cardiovascular disease in chronic dialysis patients. Nephrol Dial Transplant. 2008
- 14. Pizzarelli F, et al. Citrate high volume on-line hemodiafiltration modulates serum Interleukin-6 and Klotho levels: the multicenter randomized controlled study "Hephaestus". J Nephrol. 2021
- 15. Paim BA, et al. Oxidative stress in hypercholesterolemic LDL (low-density lipoprotein) receptor knockout mice is associated with low content of mitochondrial NADP-linked substrates and is partially reversed by citrate replacement. Free Radic Biol Med. 2008
- 16. Flythe JE, et al. Association of mortality risk with various definitions of intradialytic hypotension. J Am Soc Nephrol. 2015
- 17. Shoji T, et al. Hemodialysis-associated hypotension as an independent risk factor for two-year mortality in hemodialysis patients. Kidney Int. 2004
- 18. Belardi JA, et al. Anticoagulation with heparin The art of simplicity even in patients on dialysis. Catheter Cardiovasc Interv. 2017
- 19. Fischer KG. Essentials of anticoagulation in hemodialysis. Hemodial Int. 2007
- 20. Davenport A. What are the options for anticoagulation needs in dialysis for patients with heparin-induced thrombocytopenia? Semin Dial. 2011
- 21. David I. Ortiz-Melo ECK. Methods of Hemodialysis Anticoagulation, Handbook of dialysis therapy 2017.
- 22. Sands JJ, et al. Effects of citrate acid concentrate (citrasate™) on heparin N requirements and hemodialysis adequacy: a multicenter, prospective noninferiority trial. Blood Purif. 2012
- 23. Vandenbosch I, et al. Strategies for asymmetrical triacetate dialyser heparin-free effective haemodialysis: the SAFE study. Clin Kidney J. 2021
- 24. Kossmann, Callan, Ahmad. Fifty-five Percent Heparin Reduction is Safe with Citrate Dialysate in Chronic Dialysis Patients. ASN's 39th Annual Renal Week Meeting
- 25. Roland E. Winkler PA, Wolfgang Paetow, Grit Waitz and Hartmut Wolf. Reduction of Heparin and Oxidative Potential by Means of Citrasate™ in High-Flux Dialysis (HFD) and Online Hemodiafiltration (olHDF) in Pre and Postdilution. https://wwwintechopencom/chapters/43095. 2012
- 26. Kossmann RJ, et al. Increased efficiency of hemodialysis with citrate dialysate: a prospective controlled study. Clin J Am Soc Nephrol. 2009
- 27. Molina Nuñez M, et al. Citrate versus acetate-based dialysate in on-line haemodiafiltration. A prospective cross-over study. Blood Purif. 2015
- 28. Portolés J, et al. Anemia in Chronic Kidney Disease: From Pathophysiology and Current Treatments, to Future Agents. Front Med (Lausanne). 2021
- 29. Ma J, et al. Associations of serum and dietary magnesium with cardiovascular disease, hypertension, diabetes, insulin, and carotid arterial wall thickness: the ARIC study. Atherosclerosis Risk in Communities Study. J Clin Epidemiol. 1995
- 30. Gussak I, et al. Sudden cardiac death in nephrology: focus on acquired long QT syndrome. Nephrol Dial Transplant. 2007
- 31. Del Giorno R, et al. Consequences of Supraphysiological Dialysate Magnesium on Arterial Stiffness, Hemodynamic Profile, and Endothelial Function in Hemodialysis: A Randomized Crossover Study Followed by a Non-Controlled Follow-Up Phase. Adv Ther. 2020
- 32. Bressendorff I, et al. The Effect of Increasing Dialysate Magnesium on Serum Calcification Propensity in Subjects with End Stage Kidney Disease: A Randomized, Controlled Clinical Trial. Clin J Am Soc Nephrol. 2018
- Sakaguchi Y, et al. Hypomagnesemia is a significant predictor of cardiovascular and non-cardiovascular mortality in patients undergoing hemodialysis. Kidney Int. 2014
- 34. He K, et al. Magnesium intake and incidence of metabolic syndrome among young adults. Circulation. 2006
- 35. Okazaki H, et al. Significant positive relationship between serum magnesium and muscle quality in maintenance hemodialysis patients. Magnes Res. 2013
- 36. Yu L, et al. Association between Serum Magnesium and Erythropoietin Responsiveness in Hemodialysis Patients: A Cross-Sectional Study. Kidney Blood Press Res. 2019
- 37. Vida C, et al. A high magnesium concentration in citrate dialysate prevents oxidative stress and damage in human monocytes in vitro. Clin Kidney J. 2021
- 38. Monardo P, et al. Updates on hemodialysis techniques with a common denominator: The personalization of the dialytic therapy. Semin Dial. 2021

Nipro Renal Care is part of Nipro Corporation Japan, a leading global healthcare company established in 1954. With over 35.000 employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharmaceutical Packaging industries.

Nipro Renal Care is a global market leader with over 6 decades providing renal solutions for dialysis and dialysis-related treatment. We specialize in developing dialysis machines, water treatment systems, and a comprehensive portfolio of disposable medical equipment.

In order to address the needs of patients, healthcare professionals, and procurement managers alike, Nipro Renal Care is driven by innovation and patient safety to offer the highest quality products that optimize time, effort, and costs.

BECAUSE EVERY LIFE DESERVES AFFORDABLE CARE



www.nipro-group.com/en-en/our-company/our-locations



