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Important: Store at a temperature no lower than +4 °C. Pyrogen-free product.
Intended purpose: Preparation of dialysis fluid used in hemodialysis and related therapies **Indication:** Treatment of patients with acute or chronic renal failure when hemodialysis is prescribed by a physician. **Clinical benefits:** Facilitate the clearance of blood from uremic toxins that are able to pass through the dialysis membrane and balance out the electrolyte disturbance.
Contraindications: There are no special contraindications for the use of NiproSol. Patients who are hypersensitive to any components used in the NiproSol should not be dialyzed with this product. **Side-effects:** There are no side-effects reported for NiproSol. **Instructions:** Acid concentrate with dilution 1:44 to be used in combination with powder bicarbonate concentrate and hemodialysis machines. To be connected to the acid concentrate container, as per the machine manufacturer's instructions. Dilution is based on the following proportion. See the product label. B is the bicarbonate concentrate and H2O is deionized water, complying with: EP ISO 23500-5, and applicable to local quality requirements.
Warnings: To be used by healthcare professionals. Home dialysis is excluded. To be used under the physician's order. Ensure that the composition corresponds to the physician's prescription. Do not use if the package is damaged. The concentrated solution must be diluted immediately before use. May be used for several treatments. Once opened, the product must be used within the following 24 hours. Close the container after every use. See table below for an indication of expected dialysis fluid volume obtained from each container volume. Do not use after the expiry date stated on the package. Any serious incident that has occurred in relation to the product should be reported to the manufacturer and the competent authority of the Member State in which the user/patient is based. Report of the Member State to the EMA and national regulatory authorities.

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