NIPROTANK-D CITRASATE™

323-C

45X

Ca²⁺ 1.5 mmol/L (3.0 mEq/L)

K⁺ 2 mmol/L (2 mEq/L)

Glucose 1 g/L

EN Acid concentrate solution for hemodialysis and related therapies

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C₆H₈O₇: 6.92 g/L

H₂O Ph. Eur. q.s.: 1000 mL

EN Electrolytic composition of the ready-to-use solution

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EN: Instructions for use

Important: Store at a temperature no lower than 4°C. Pyrogen-free product
Intended purpose: Preparation of dialysis fluid used in hemodialysis and related therapies. **Indications:** Treatment of patients with acute or chronic renal failure when hemodialysis is prescribed by a physician. **Clinical Benefits:** Facilitate the clearance of blood from uremic toxins that are able to pass through the dialysis membrane and balance out the electrolytic disturbance. **Contraindications:** There are no special contraindications for the use of Niprotank-D Citrasate. Patients who are hypersensitive to any components used in the Niprotank-D Citrasate should not be dialyzed with this product. **Side-effects:** There are no side-effects reported for Niprotank-D Citrasate. **Instructions:** Acid concentrate with dilution 1:44 to be used in combination with powder bicarbonate concentrate and hemodialysis machines. To be connected to the acid concentrate connector of the central concentrate supply system. The acid concentrate from the central concentrate supply system must be used with a dialysis machine, following the instructions of the machine manufacturer. Dilution is based on the following properties. See the product label. B is the bicarbonate concentrate and H2O is deionized water, complying with EP ISO 23500-5, and applicable local quality requirements. **Warnings:** To be used by healthcare professionals. Home dialysis is excluded. To be used under the physician's order. Ensure that the composition corresponds to the physician's prescription. Do not use if package is damaged. The concentrated solution must be diluted immediately before use. Once opened, the product must be used continuously without disconnecting from the central concentrate supply system. Do not use after disconnection from the central concentrate supply system. Do not use after the expiry date stated on the package. Citrate, a component of this product, may reduce the need for heparin and may lower plasma Ca and Mg concentrations. See table below for an indication of expected dialysis fluid volume obtained from each container volume. Any serious incident that has occurred in relation to the product should be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is based. Dispose of the product according to local and national regulations.

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Theramed AG Therapie und Medizintechnik
UK Responsible Person: Nipro Medical UK Ltd.