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☞ Instructions for use Acid Concentrate and Citrasate® for Bicarbonate Hemodialysis

Indications (ready-to-use solution): Acute and chronic renal failure, hyperhydration, intoxication, correction of acid/base metabolism and electrolyte status, adjustment of blood/plasma/body temperature. The clinical data reviewed suggest that Citrate/Citrasate® should be favorable for the following types of patient indications: hemodynamic instability, low PTH or a dynamic bone disease; hypercalcemia; vascular calcifications; patients with an indication to reduced heparin or heparin-free dialysis; better correction of strong metabolic acidosis; hypermagnesiemia; difficulty in reaching the prescribed dose of dialysis, malnutrition, reduction of inflammation and oxidative stress. **Type of application:** For extracorporeal bicarbonate dialysis. Use acid concentrate only in combination with bicarbonate hemodialysis concentrate (NaHCO3 8,4% or 6,6%) or the bicarbonate cartridge in the given dilution; including on-line procedures, but only in combination with a hemodialysis monitor capable and intended for online treatments (see relevant dialysis machine operator manual). Perform a gas analysis on the ready-to-use bicarbonate hemodialysis solution prior to use, employing an appropriate method. **Side effects:** Hypotonia, hypertension, arrhythmias, nausea, vomiting, spasms, chills, headache dizziness, numbness, palpitations, fatigue, malaise, itching have been reported. **Other side effects:** Alkalosis resulting in cardiac arrhythmia, paresthesia, cramps; reduced cerebral blood flow; respiratory distress, headache and a procalcifying effect; inflammation resulting in no calcium homeostasis; inadequate magnesium or calcium plasma concentration resulting in cramps, arrhythmias and hypotension, cardiac complications; vascular calcifications and change of blood pressure and phosphate concentrations. Bleeding or clotting during treatment may result in low dialysis efficacy, hemorrhage and vascular access problems. **Contraindications:** Hyperkalemia (applicable for potassium containing concentrates), hypokalemia (applicable for potassium-free concentrates), and unmanageable coagulation disorders. For patients with strong hemodynamic instability an alternative method of extracorporeal treatment may be indicated based on the overall patient indications. **Notice:** Follow the manufacturer instructions for the dialysis machine. The concentration of the ready-made dialysis solution must be precisely monitored. Check the conductivity and the pH level of the dialysate before dialysis treatment and each time a new concentrate is used with the machine. Acid and basic concentrate for bicarbonate hemodialysis should only be diluted with water of a suitable quality according to ISO 23500-3 immediately before use. Keep the container closed when not in use. Residues should be discarded. Use only undamaged containers. Depending on the volume, the concentrate can be used for more than one dialysis treatment. Recommendation for use of opened cases: canister max. 1 week, container/bag in box max. 3 months. Do not use after the expiration date. Do not use already opened containers after the expiration date. The concentrate is endotoxin-free. For safe application and disposal, see safety data sheet of manufacturer. Check-up of metabolizing capacity of citrate for all patients. Patients with small muscle mass and/or restricted liver function may not be suited without further correction of concentration of free calcium concentration. Compared with acetate dialysate, notably higher calcium and magnesium concentrations, may be required in the citrate-containing dialysate to achieve similar systemic post-dialytic values. Because one citrate ion will be converted quickly during dialysis treatment into 3 bicarbonate ions, a change of acid/base balance takes place with the possible result of a temporary alkalosis. Therefore, the acid/base balance during treatment has to be monitored and corrected by changing the pre-dialysis bicarbonate concentration or bicarbonate profiling during treatment. During use of Citrasate® without reduction of anticoagulants may pose a risk for bleeding and hypercoagulation, which in turn may be associated with increased venous access complications. Any serious incident that has occurred in connection with the product must be notified to the manufacturer and the competent authority of the Member State in which the user is established. For safe connection of the containers or bag-in-box systems to central ring line systems, the instructions in the appendix to the instructions for use must be followed!

UK Responsible Person:

Nipro Medical UK, Ltd Units 12-14 South Point, Ensign Way, Hamble-Le-Rice, Southampton, Hampshire, SO31 4RF, United Kingdom

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