



0800 08 588 08

Call Nipro Customer Service for more information on products or for help in switching patients to the 4SURE portfolio.

Nipro Diabetes Care is part of Nipro Corporation Japan, a leading global healthcare company established in 1954. With over 28,000 employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharmaceutical Packaging industries.

Nipro Diabetes Care is a global supplier of high quality products for people living with type 1 or type 2 diabetes.

Our portfolio of meters, test strips, pen needles, and lancets capture the essentials of self-testing in the most simple, accurate, and accessible form. This way we ensure the best life possible!

BECAUSE EVERY LIFE DESERVES AFFORDABLE CARE



Nipro Diagnostics UK Ltd. :
Units 12-14 South Point, Ensign Way, Hamble, Hampshire SO31 4RF, United Kingdom.
diagnostics-uk@nipro-group.com | www.nipro-diagnostics.co.uk | 0800 08 588 08

**LET US HELP YOU
MANAGE DIABETES**



DIABETES CARE

Bro-4Sure_PenNeedles - EN-UK - 12.Mar2019



PEN NEEDLES
ADVANCED FEATURES FOR COMFORTABLE USE



SURE™ PEN NEEDLES

Sure and simple

ADVANCED FEATURES

Thin wall technology

Large inner diameter allows for optimum insulin flow. This translates to reduced injection force, reduced injection duration, and - most importantly - reduced pain.

Double silicone coating

Special double layer silicone coating ensures a smooth needle insertion, resulting in minimal pain.

Ultra sharp tip

Ultra sharp tip ensures a clean, precise penetration of the skin.

COMFORTABLE USE

Clear needle cap

Clear cap to minimise needle stick injuries and to aid in air shot inspection.

Universal fit

Compatible with all major insulin pens.

Excellent ergonomic grip

For easy cap removal.



YOUR
COST-EFFECTIVE
CHOICE



SWITCH TO SURE™ TODAY

A wide portfolio for an easy switch

Reference code	Length	Diameter	Color
4S-IPN-0671-31-4	4 mm	31 G (0.25 mm)	Light Green
4S-IPN-0671-32-4	4 mm	32 G (0.23 mm)	Medium Green
4S-IPN-0671-31-5	5 mm	31 G (0.25 mm)	Purple
4S-IPN-0671-32-5	5 mm	32 G (0.23 mm)	Dark Purple
4S-IPN-0671-31-6	6 mm	31 G (0.25 mm)	Light Blue
4S-IPN-0671-32-6	6 mm	32 G (0.23 mm)	Dark Blue
4S-IPN-0671-31-8	8 mm	31 G (0.25 mm)	Orange

LET US HELP YOU MANAGE DIABETES