

# ACID CONCENTRATES

BAGS, CANISTERS, OCTABINS



# ACID CONCENTRATES

BAGS, CANISTERS, OCTABINS

Ready-to-use liquid acid concentrates intended for use with powder bicarbonate concentrate. Acid concentrate is an essential component in the preparation of dialysis fluid, whereby the hemodialysis machine mixes liquid acid concentrate with powder bicarbonate concentrate and purified water.

The standard mixing ratio of Nipro acid concentrates is 1:44.



## PORTFOLIO

	Na <sup>+</sup>	Ca <sup>++</sup>	Mg <sup>++</sup>	K <sup>+</sup>	Cl <sup>-</sup>	CH <sub>3</sub> COO <sup>-</sup>	HCO <sub>3</sub> <sup>-</sup>	Glucose	Osmolarity
	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	g/l	mOsm/l
306-A (A+B)	139	1,25	0,48	1,50	103,14	2,85	37,96	1,50	294,50
307-A (A+B)	139	1,75	0,48	1,50	104,14	2,85	37,96	1,50	296,00
308-A (A+B)	139	1,50	0,48	1,50	103,64	2,85	37,96	1,00	292,48
310-A (A+B)	139	1,50	0,50	1,00	103,20	2,85	37,96	1,00	291,56
313-A (A+B)	140	1,50	0,50	1,50	108,50	3,00	34,00	1,00	294,55
314-A (A+B)	140	1,25	0,50	1,50	108,00	3,00	34,00	1,00	293,80
316-A (A+B)	140	0,63	0,25	1,50	102,44	4,00	36,81	1,00	291,18
318-A (A+B)	140	1,75	0,50	2,00	109,50	3,00	34,00	1,50	299,08
319-A (A+B)	139	1,75	0,48	3,00	105,64	2,85	37,96	1,50	299,00
320-A (A+B)	139	1,25	0,50	3,00	104,70	2,85	37,96	1,50	297,59
321-A (A+B)	140	1,75	0,50	2,00	109,50	3,00	34,00	1,00	296,30
323-A (A+B)	140	1,50	0,50	2,00	109,00	3,00	34,00	1,00	295,55
324-A (A+B)	140	1,25	0,50	2,00	108,50	3,00	34,00	1,00	294,80
326-A (A+B)	139	1,50	0,48	2,50	104,64	2,85	37,96	1,50	297,25
327-A (A+B)	140	1,50	0,50	3,00	110,00	3,00	34,00	1,00	297,55
330-A (A+B)	140	1,75	0,50	3,00	110,50	3,00	34,00	1,00	298,30
331-A (A+B)	140	1,25	0,50	3,00	109,50	3,00	34,00	1,00	296,80
332-A (A+B)	140	1,50	0,50	4,00	111,00	3,00	34,00	1,00	299,55

## PACKAGING DETAILS

	Volume (L)	Reference code	Units per box	Units per pallet
Bags	3,8	BSOL45-*A4LBX	3	144
	4,7	BSOL45-*A5LBX	2	128
Stackable canister(s)	3,9	SOL45-*A4L	N/A	132
	5	SOL45-*A5L	N/A	132
	10,5	SOL45-*A11L	N/A	60
Canister(s) in a box	3,9	SOL45-*A4LBX	2	128
	5	SOL45-*A5LBX	2	128
Octabin	500	CSOL45-*A500L	N/A	1
	600	CSOL45-*A600L	N/A	1
Big box stackable canisters	3,9	SOL45-*A4LBB	96	96
	5	SOL45-*A5LBB	96	96

\* acid concentrate formula

Not all products are registered and approved for sale in all countries or regions. Indications of use may also vary by country and region. Please contact your country representative for product availability and information.


**SHELF LIFE:** 2 years

Nipro Renal Care is part of Nipro Corporation Japan, a leading global healthcare company established in 1954. With over 28.000 employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharmaceutical Packaging industries.

Nipro Renal Care is a global market leader with over 5 decades providing renal solutions for dialysis and dialysis-related treatment. We specialize in developing dialysis machines, water treatment systems, and a comprehensive portfolio of disposable medical equipment.

In order to address the needs of patients, healthcare professionals, and procurement managers alike, Nipro Renal Care is driven by innovation and patient safety to offer the highest quality products that optimize time, effort, and costs.

**BECAUSE EVERY LIFE DESERVES AFFORDABLE CARE**



**Nipro Medical Europe** : European Headquarters, Blokhuisstraat 42, 2800 Mechelen, Belgium  
T: +32 (0)15 263 500 | F: +32 (0)15 263 510 | [medical@nipro-group.com](mailto:medical@nipro-group.com) | [www.nipro-group.com](http://www.nipro-group.com)

Please contact your local representative for more information.

**NIPRO MEDICAL AUSTRIA GMBH :**

Divischgasse 4, 1210 Wien, AUSTRIA | T: +43 1 532 23 14 | F: +43 1 532 23 14 89

**NIPRO EUROPE - EGYPT :**

Nile City Towers, 22nd Floor, North Tower, Nile City, Towers, Cornich El Nile, 11624 Ramelt Beaulac, Cairo, EGYPT

**NIPRO MEDICAL FRANCE S.A. :**

Biopôle Clermont-Limagne, 63360 Saint Beauzire, FRANCE | T: +33 (0)473 33 41 00 | F: +33 (0)473 33 41 09

**NIPRO MEDICAL GERMANY GMBH :**

Kokkolastrasse 5, 40882 Ratingen, GERMANY | T: +49 (0)2102 564 60 81 | F: +49 (0)2102 564 60 90

**NIPRO MEDICAL EUROPE NV - ITALY :**

Centro Direzionale Milanofiori, Strada 1 - Palazzo F1, 20090 Assago (Milano), ITALY | T: +39 (0)2 57 50 00 57 | F: +39 (0)2 57 51 81 11

**NIPRO MEDICAL IVORY COAST :**

Selmer Lot n°18, Abidjan-Cocody Riviera, Abidjan, Ivory Coast

**NIPRO MEDICAL EUROPE NV - NETHERLANDS :**

Verlengde Poolseweg 16, 4818 CL Breda, NETHERLANDS | T: +31 (0)76 524 50 99 | F: +31 (0)76 524 46 66

**NIPRO MEDICAL NIGERIA LTD. :**

9 Adelabu Close, Off Toyin Street, 100271 Ikeja, Lagos State, NIGERIA | T: +234 (0)802 706 7065

**NIPRO MEDICAL POLAND SP. Z.O.O. :**

Ul. Panska 73, 00-834 Warszawa, POLAND | T: +48 (0)22 31 47 155 | F: +48 (0)22 31 47 152

**NIPRO EUROPE - PORTUGAL :**

Avd. Da Liberdade 249, 1º Andar, 1205-143 Lisboa, PORTUGAL | T: +34 (0)91 878 29 21 | F: +34 (0)91 878 28 40

**NIPRO MEDICAL EUROPE - RUSSIAN REPRESENTATIVE OFFICE :**

12 Krasnopresnenskaya Nab., Office 1407, entrance 6, 123610 Moscow, RUSSIA | T: +7 (0)495 258 1364 | F: +7 (0)495 258 1365

**NIPRO MEDICAL SENEGAL S.U.A.R.L. :**

27 Avenue Georges Pompidou, Dakar, SENEGAL

**NIPRO MEDICAL D.O.O. BEOGRAD :**

Bastovanska 68a, 11000 Belgrade, SERBIA | T: +381 (0)11 75 15 578

**NIPRO MEDICAL SOUTH AFRICA (PTY) LTD :**

4B Dwiyka Street, Stikland Industria, Cape Town, 7530, SOUTH AFRICA | T: +27 21 949-2635 | F: +27 21 949-2397

Unit 20&21, Falcon Lane, Lanseria Business Park, Erf 805 Lanseria Corporate Estate, Pelindaba Rd.,

Lanseria Ext. 26, Gauteng, SOUTH AFRICA | T: +27 11 431 1114 / 26 | F: +27 11 431 1115

**NIPRO MEDICAL EUROPE NV - SPAIN :**

Potigono Los Frailes n° 93 y 94, Daganzo, 28814 Madrid, SPAIN | T: +34 (0)91 884 5531 | F: +34 (0)91 878 2840

**NIPRO MEDIKAL SAĞ.HIZ.TIC.LTD.ŞTİ. :**

Aziyiyeh Mah.Pilot Sok.17/A, 6540 Çankaya-Ankara, TURKEY | T: +90 (0)312 442 21 12 | F: +90 (0)312 442 21 92

**NIPRO MEDICAL UK LTD. :**

25 Barnes Wallis Road, Segensworth East, Fareham Hampshire PO15 5TT, UNITED KINGDOM | T: +44 148 985 48 30