



3-PART SYRINGE WITH NEEDLE

STERILE, SINGLE-USE SYRINGE WITH DETACHABLE NEEDLE
FOR INJECTION OR ASPIRATION



3-PART SYRINGE WITH NEEDLE

STERILE, SINGLE-USE SYRINGE WITH DETACHABLE NEEDLE FOR INJECTION OR ASPIRATION

- Ultra sharp three-beveled hypodermic needle
- Transparent barrel for easy recognition of air bubbles
- Smooth plunger motion along the entire length of the barrel
- Large finger grips for secure handling
- Latex-free, PVC-free
- Available in centric and eccentric tip
- Available in 18-27 G
- For use by healthcare professionals only

CE₀₁₂₃

Class II Medical Devices
Rule 6
MDD 93/42/EEC
UMDNS 13929

PRODUCT COMPLIANCE

- CE-marked, Class II Medical Devices, Rule 6, MDD 93/42/EEC, UMDNS 13929
- Complies with the following norms, directives, and regulations:
 - ISO 594-1
 - ISO 7864
 - ISO 7886-1
 - ISO 9626
 - EN ISO 13485
 - EN ISO 15223-1
 - ISO 11607-1
- Labels contain 15 languages

MANUFACTURING DETAILS

Legal manufacturer: Nipro India Corporation PVT. LTD.

Country of origin: India

STERILIZATION AND SHELF LIFE

Sterilization

Gamma-ray



Single-use only



Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a yellow sticker that turns red when sterilized.

Shelf life

5 years



MATERIALS USED



Latex-free
PVC-free

PRODUCT RANGE OVERVIEW

Reference code	Description	Size (ml)	Type	Scaling (ml)	Dead Volume (ml)
SY3-1SCTU-XXXX-GIN2	3 part syringe 1 ml luer slip	1	Centric	0,01	0,07
SY3-2SC-XXXX-GIN2	3 part syringe 2,5 ml luer lock	2,5	Centric	0,1	0,07
SY3-3SC-XXXX-GIN2	3 part syringe 3 ml luer lock	3	Centric	0,1	0,07
SY3-5SC-XXXX-GIN2	3 part syringe 5 ml luer slip	5	Centric	0,2	0,07
SY3-10SC-XXXX-GIN2	3 part syringe 10 ml luer slip	10	Centric	0,2	0,10
SY3-20ESC-XXXX-GIN2	3 part syringe 20 ml luer slip	20	Eccentric	1	0,15
SY3-20LC-XXXX-GIN2	3 part syringe 20 ml luer lock	20	Centric	1	0,15
SY3-50ESC-XXXX-GIN2	3 part syringe 50 ml luer slip	50	Eccentric	1	0,20
SY3-50LC-XXXX-GIN2	3 part syringe 50 ml luer lock	50	Centric	1	0,20

XXXX: The first two digits express the needle diameter (G), and the last two digits express the needle length (mm).

PACKAGING DETAILS

Reference	Size (ml)	Scaling	Inner box					Outer box				
			Quantity	L (mm)	W (mm)	H (mm)	Weight (kg)	Quantity	L (cm)	W (cm)	H (cm)	Weight (kg)
SY3-1SCTU-XXXX-GIN2	1	0,01	100	340	155	75	0,51	1000	385	370	350	6,2
SY3-2SC-XXXX-GIN2	2,5	0,1	100	340	145	90	0,52	1000	485	330	350	6,5
SY3-3SC-XXXX-GIN2	3	0,1	100	340	145	90	0,55	1000	485	330	350	6,7
SY3-5SC-XXXX-GIN2	5	0,2	100	425	155	115	0,83	1500	585	545	435	14,9
SY3-10SC-XXXX-GIN2	10	0,2	100	410	185	150	1,18	800	590	400	460	11,3
SY3-20ESC-XXXX-GIN2	20	1	50	415	175	115	0,95	500	585	440	425	11,6
SY3-20LC-XXXX-GIN2	20	1	50	415	175	115	-	500	585	440	425	-
SY3-50ESC-XXXX-GIN2	50	1	30	400	183	135	-	240	550	380	410	-
SY3-50LC-XXXX-GIN2	50	1	30	400	183	135	-	240	550	380	410	-

Languages

15 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), Finnish (FI), Polish (PL), Hebrew (HE), and Arabic (AR)


Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE R	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	✓	✓	✓
	—	—	✓

