



SYRINGE WITH NEEDLE

STERILE, SINGLE-USE SYRINGE WITH DETACHABLE
NEEDLE FOR INJECTION



SYRINGE WITH NEEDLE

STERILE, SINGLE-USE SYRINGE WITH DETACHABLE NEEDLE FOR INJECTION

- Transparent barrel for precise volume reading
- Smooth plunger motion along entire length of barrel
- Large finger grips for secure handling
- Latex-free, BPA-free, DEHP-free, PVC-free
- Available in centric or eccentric tip
- Available with luer lock or luer slip
- Available in 14-30 G
- For use by healthcare professionals only.



Class IIa Medical Devices
Rule 6
MDD 93/42/EEC
UMDNS 13940

PRODUCT COMPLIANCE

- CE-marked, Class IIa Medical Devices, Rule 6, MDD 93/42/EEC, UMDNS 13940
- Complies with the following norms, directives, and regulations:
 - ISO 594-1
 - EN 1707
 - ISO 7864
 - ISO 7886-1
 - ISO 13485
 - ISO 10993-1
 - ISO 10993-4
 - ISO 10993-4
 - ISO 10993-5
 - ISO 10993-7
 - ISO 11607-1
 - ISO 11607-2
 - EN ISO 13485
 - ISO 15223-1
- Labels contain 15 languages

MANUFACTURING DETAILS

Legal manufacturer: Shanghai Kindly Enterprise Development Group Co., Ltd.

Country of origin: China

STERILIZATION AND SHELF LIFE

Sterilization

EtO (Ethylene oxide)



Single-use only




Shelf life

5 years



MATERIALS USED



DEHP-free
BPA-free
Latex-free
PVC-free

PRODUCT RANGE OVERVIEW

Reference code	Description	Size (ml)	Tip type
SY3-1SC-XXXX-EC	3 part syringe 1 ml luer slip	1	CENTRIC
SY3-2SC-XXXX-EC	3 part syringe 2 ml luer slip	2	CENTRIC
SY3-2LC-XXXX-EC	3 part syringe 2 ml luer lock	2	CENTRIC
SY3-3SC-XXXX-EC	3 part syringe 3 ml luer slip	3	CENTRIC
SY3-5SC-XXXX-EC	3 part syringe 5 ml luer slip	5	CENTRIC
SY3-5LC-XXXX-EC	3 part syringe 5 ml luer lock	5	CENTRIC
SY3-10SC-XXXX-EC	3 part syringe 10 ml luer slip	10	CENTRIC
SY3-10ESC-XXX-EC	3 part syringe 10 ml luer slip	10	ECCENTRIC
SY3-10LC-XXXX-EC	3 part syringe 10 ml luer lock	10	CENTRIC
SY3-20SC-XXXX-EC	3 part syringe 20 ml luer slip	20	ECCENTRIC
SY3-20LC-XXXX-EC	3 part syringe 20 ml luer lock	20	CENTRIC
SY3-30SC-XXXX-EC	3 part syringe 30 ml luer slip	30	ECCENTRIC
SY3-30LC-XXXX-EC	3 part syringe 30 ml luer lock	30	CENTRIC

XXXX: The first two digits express the needle diameter (G), and the last two digits express the needle length (mm).

PACKAGING DETAILS

Reference	Size (ml)	Inner box				Outer box				Net weight (kg)	Gross weight (kg)
		Quantity	L (mm)	W (mm)	H (mm)	Quantity	L (cm)	W (cm)	H (cm)		
SY3-1SC-XXXX-EC	1	100	210	140	108	3200	44	46	58	12,00	16,6
SY3-2SC-XXXX-EC	2	100	220	207	90	2400	58	44	46	11,5	15,70
SY3-2LC-XXXX-EC	2	100	220	207	90	2400	58	44	46	11,5	15,70
SY3-3SC-XXXX-EC	3	100	220	207	90	2400	58	44	46	11,5	15,70
SY3-5SC-XXXX-EC	5	100	420	186	72	1800	58	44	46	13,2	18,50
SY3-5LC-XXXX-EC	5	100	420	186	72	1800	58	44	46	13,2	18,50
SY3-10SC-XXXX-EC	10	100	420	215	112	1000	58	44	46	10,10	13,60
SY3-10ESC-XXX-EC	10	100	420	215	112	1000	58	44	46	13,70	17,80
SY3-10LC-XXXX-EC	10	100	420	215	112	1000	58	44	46	12,50	16,70
SY3-20SC-XXXX-EC	20	50	215	210	140	800	58	44	46	12,70	16,40
SY3-20LC-XXXX-EC	20	50	215	210	140	800	58	44	46	12,70	16,40
SY3-30SC-XXXX-EC	30	50	285	210	153	600	58	44	46	12,20	15,80
SY3-30LC-XXXX-EC	30	50	285	210	153	600	58	44	46	12,20	15,80

Languages

15 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), Finnish (FI), Russian (RU), Turkish (TR), and Polish (PL)

Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
REF	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
	—	—	✓
	✓	✓	✓

