




## STANDARD SCALP VEIN SET

FOR INFUSION THERAPY



# STANDARD SCALP VEIN SET

FOR INFUSION THERAPY

- Ultra sharp three-beveled needle
- Siliconization of needle minimizes penetration and gliding force
- Large soft wings for easy grip
- Color-coded to conform to ISO 6009
- Latex-free, DEHP-free
- For use by healthcare professionals only

CE

Class IIa Medical Devices  
Rule 7  
MDD 93/42/EEC  
UMDNS 12752

## PRODUCT COMPLIANCE

- CE-marked, Class IIa Medical Devices, Rule 7, MDD 93/42/ECC, UMDNS 12752
- Complies with the following norms, directives, and regulations:
  - ISO 594-1
  - ISO 594-2
  - CSN EN 1707
  - ISO 6009
  - ISO 7864
  - ISO 11607-1
  - ISO 11607-2
  - ISO 13485
  - ISO 15223-1
- Labels contain 12 languages

## MANUFACTURING DETAILS

**Legal manufacturer:** Nipro Corporation

**Country of origin:** Thailand

## STERILIZATION AND SHELF LIFE

### Sterilization

EtO (Ethylene oxide)



### Single-use only



### Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that ensures sterility. This indicator comes in the form of a blue sticker that turns red when sterilized.

### Shelf life

5 years



# MATERIALS USED



**DEHP-free**  
**Latex-free**

# PRODUCT RANGE OVERVIEW

Reference code	Description	Tubing length	Purpose	Color	Flow rate not less than	Priming volume
PSV-18-ET	Winged Needle Set 18Gx3/4" (1.2x19 mm)	300 mm	IV/SC	Pink	46,5 ml/min	0,46 ml
PSV-19-ET	Winged Needle Set 19Gx3/4" (1.1x19 mm)	300 mm	IV/SC	Orange	36,5 ml/min	0,47 ml
PSV-20-ET	Winged Needle Set 20Gx3/4" (0.9x19 mm)	300 mm	IV/SC	Yellow	25,8 ml/min	0,46 ml
PSV-21-ET	Winged Needle Set 21Gx3/4" (0.8x19 mm)	300 mm	IV/SC	Green	20 ml/min	0,45 ml
PSV-22-ET	Winged Needle Set 22Gx3/4" (0.7x19 mm)	300 mm	IV/SC	Black	12 ml/min	0,44 ml
PSV-23-ET	Winged Needle Set 23Gx3/4" (0.6x19 mm)	300 mm	IV/SC	Blue	7,8 ml/min	0,44 ml
PSV-24-ET	Winged Needle Set 23Gx3/4" (0.6x19 mm)	300 mm	IV/SC	Purple	5,6 ml/min	0,44 ml
PSV-25-ET	Winged Needle Set 25Gx3/4" (0.5x19 mm)	300 mm	IV/SC	Orange	3,7 ml/min	0,45 ml
PSV-27-ET	Winged Needle Set 27Gx3/4" (0.4x19 mm)	300 mm	IV/SC	Grey	0,35 ml/min	0,45 ml

IV: Intravenous, SC: Subcutaneous, E: Ethylene Oxide, T: Thailand

# PACKAGING DETAILS

Box	Quantity (pcs)	Net weight (kg)	Gross weight (kg)	Length (mm)	Width (mm)	Height (mm)	m <sup>3</sup>
Inner box	50	—	—	247	110	78	—
Outer box	1000	5,40	6,40	572	263	344	0,052

## Languages

12 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), and Finnish (FI)







## Transport conditions

Closed and dry

## Storage conditions

Open the packaging only immediately before use to guarantee sterility. Store in a dry place away from direct heat and sunlight.

## Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	—	✓	✓
<b>REF</b>	—	✓	✓
<b>LOT</b>	✓	✓	✓
	✓	✓	✓
<b>STERILE</b>	✓	✓	✓
<b>CE</b>	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
<b>EC REP</b>	✓	✓	✓
	—	✓	✓

