



3-PART SYRINGE WITHOUT NEEDLE

STERILE, SINGLE-USE SYRINGE FOR INJECTION OR ASPIRATION



3-PART SYRINGE WITHOUT NEEDLE

STERILE, SINGLE-USE SYRINGE FOR INJECTION OR ASPIRATION

- Transparent barrel for easy recognition of air bubbles
- Smooth plunger motion along the entire length of the barrel
- Large finger grips for secure handling
- Not made with latex or PVC
- Available in centric or eccentric tip
- Available in luer lock or luer slip

CE₀₁₂₃

Class I Medical Devices
Rule 2
MDD 93/42/EEC
UMDNS 13929

PRODUCT COMPLIANCE

- CE-marked, Class I Medical Devices, Rule 2, MDD 93/42/EEC, UMDNS 13929
- Complies with the following norms, directives, and regulations:
 - ISO 594-1
 - ISO 7864
 - ISO 7886-1
 - ISO 9626
 - ISO 11607-1
 - EN ISO 13485
 - EN ISO 15223-1
- Labels contain 15 languages

MANUFACTURING DETAILS

Legal manufacturer: Nipro India Corporation PVT. LTD.

Country of origin: India

STERILIZATION AND SHELF LIFE

Sterilization

Gamma-ray



Single-use only



Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a yellow sticker that turns red when sterilized.

Shelf life

5 years



MATERIALS USED



Not made with latex or PVC

PRODUCT RANGE OVERVIEW

Reference code	Description	Size (ml)	Type	Scaling (ml)	Dead Volume (ml)
SY3-1SCTU-GIN2	3 part syringe 1ml luer slip	1	Centric	0,01	0,07
SY3-2SC-GIN2	3 part syringe 2.5 ml luer lock	2,5	Centric	0,1	0,07
SY3-5SC-GIN2	3 part syringe 5ml luer slip	5	Centric	0,2	0,07
SY3-10SC-GIN2	3 part syringe 10ml luer slip	10	Centric	0,2	0,10
SY3-20ESC-GIN2	3 part syringe 20 ml luer slip	20	Eccentric	1	0,15
SY3-20LC-GIN2	3 part syringe 20ml luer lock	20	Centric	1	0,15
SY3-50LC-GIN2	3 part syringe 50ml luer lock	50	Centric	1	0,20

PACKAGING DETAILS

Reference	Size (ml)	Inner box					Outer box				
		Quantity	L (mm)	W (mm)	H (mm)	Weight (kg)	Quantity	L (cm)	W (cm)	H (cm)	Weight (kg)
SY3-1SCTU-GIN2	1	100	400	135	65	0,40	2000	550	335	410	8,9
SY3-2SC-GIN2	2.5	100	340	120	95	0,40	2000	500	490	360	9,8
SY3-5SC-GIN2	5	100	425	125	115	0,68	1500	590	390	440	12,1
SY3-10SC-GIN2	10	100	450	145	145	1,00	1200	600	445	470	14,1
SY3-20ESC-GIN2	20	50	415	175	115	0,87	500	585	440	425	10,2
SY3-20LC-GIN2	20	50	415	175	115	0,87	500	585	360	425	10,2
SY3-50LC-GIN2	50	30	400	183	135	1,05	240	550	380	410	10,0

Languages

15 languages on all levels of packaging:

English (EN), French (FR), Dutch (NL), German (DE), Spanish (ES), Italian (IT), Portuguese (PT), Greek (EL), Swedish (SV), Danish (DA), Norwegian (NO), Finnish (FI), Polish (PL), Hebrew (HE), and Arabic (AR)

Transport conditions

Closed and dry

Storage conditions

Open the packaging only immediately before use to guarantee sterility.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pieces)	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE R	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
EC REP	✓	✓	✓
	✓	✓	✓
	—	—	✓

