



HUBER NEEDLE (Straight)

STERILE, SINGLE-USE, NON-CORING HUBER NEEDLE FOR
ADMINISTRATION OF FLUIDS AND DRUGS OR BLOOD SAMPLING



HUBER NEEDLE (Straight)

STERILE, SINGLE-USE, NON-CORING HUBER NEEDLE FOR ADMINISTRATION OF FLUIDS AND DRUGS OR BLOOD SAMPLING

- To be inserted into the implanted vascular port
- Indicated for the intra-vascular administration of fluids and drugs or blood sampling through the implanted port
- Non-coring needle with cap
- Straight needle
- Not made with latex, DEHP, or PVC
- Available in 20, 22 G
- Available in 25, 32, 38 mm (1, 1 ¼, 1 ½ inches, respectively)

CE 0123

Class IIb Medical Devices

Rule 7

UMDNS 12729

GMDN 17701

Product Compliance

- Complies with the following norms, directives, and regulations:
 - ISO 7864
 - ISO 80369-7
 - ISO 10555-3
 - ISO 8536-4

Manufacturing DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

Sterilization and Shelf Life

Sterilization

EtO (Ethylene oxide)



Single-use only



Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a violet sticker that turns to blue when sterilized.

Shelf life

3 years



Materials Used



PRODUCT RANGE OVERVIEW

Reference code	Description	Needle outer diameter	Needle length
HCN-S-2025-ET	Huber coreless Ndl straight 20 G 25 mm	20 G	25 mm
HWN-S-2032-ET	Huber coreless Ndl straight 20 G 32 mm	20 G	32 mm
HCN-S-2038-ET	Huber coreless Ndl straight 20 G 38 mm	20 G	38 mm
HCN-S-2225-ET	Huber coreless Ndl straight 22 G 25 mm	22 G	25 mm
HCN-S-2238-ET	Huber coreless Ndl straight 22 G 38 mm	22 G	38 mm

Needle cap lengths:

25 mm (1") needle length: 30 mm needle cap

32 mm (1 1/4") and 38 mm (1 1/2") needle length: 42,7 mm needle cap

Packaging Details

Individual unit			Inner box				Outer box					
L (mm)	H (mm)	Weight (g)	Quantity (pcs)	L (mm)	W (mm)	H (mm)	Weight (g)	Quantity (pcs)	L (mm)	W (mm)	H (mm)	Weight (kg)
85	65	0,7	50	195	90	70	120	1000	414	387	211	3,2

Languages

26 languages of the IFU available:







(BG) Bulgarian, (CS) Czech, (DA) Danish, (DE) German, (EL) Greek, (EN) English, (ES) Spanish, (ET) Estonian, (FI) Finnish, (FR) French, (HR) Croatian, (HU) Hungarian, (IT) Italian, (LT) Lithuanian, (LV) Latvian, (NL) Dutch, (NO) Norwegian, (PL) Polish, (PT) Portuguese, (RO) Romanian, (RU) Russian, (SK) Slovak, (SL) Slovenian, (SR) Serbian, (SV) Swedish, (TR) Turkish.

Transport and storage conditions

Keep dry and away from sunlight.

Label details

	Individual unit	Inner box	Outer box
Product name	✓	✓	✓
Product size	✓	✓	✓
Quantity (pcs)	—	✓	✓
REF	✓	✓	✓
LOT	✓	✓	✓
	✓	✓	✓
STERILE EO	✓	✓	✓
CE	✓	✓	✓
	✓	✓	✓
	✓	✓	✓
	✓	✓	✓

	Individual unit	Inner box	Outer box
	✓	✓	✓
EC REP	✓	✓	✓
	—	✓	✓
	—	—	✓
	—	—	✓
	—	—	✓
	✓	✓	✓
Not made with LATEX	✓	✓	✓
Not made with DEHP	✓	✓	✓

