



HUBER NEEDLE (WITH 90° ANGLE)

STERILE, SINGLE-USE, NON-CORING HUBER NEEDLE FOR INFUSION
AND BLOOD SAMPLING



HUBER NEEDLE (with 90° angle)

STERILE, SINGLE-USE, NON-CORING HUBER NEEDLE FOR INFUSION AND BLOOD SAMPLING

- Intended for the infusion of drug solutions into or blood sampling from a reservoir implanted in the body
- Non-coring needle with cap
- Needle fixed at 90° angle
- Not made with latex, DEHP, or PVC
- Available in 19, 20, 22 G
- Available in 19, 25, 32, 38 mm (¾, 1, 1 ¼, 1 ½ inches, respectively)

CE 0123

Class IIb Medical Devices

Rule 7

UMDNS 12729

GMDN 17701

Product Compliance

- Complies with the following norms, directives, and regulations:
 - ISO 7864
 - ISO 80369-7
 - ISO 10555-3
 - ISO 8536-4

Manufacturing DETAILS

Legal manufacturer: Nipro Corporation

Country of origin: Thailand

Sterilization and Shelf Life

Sterilization

EtO (Ethylene oxide)



Single-use only



Indicator sticker


Each outer box (or shipping carton) contains a chemical indicator that indicates sterility. This indicator comes in the form of a violet sticker that turns to blue when sterilized.

Shelf life

3 years



Materials Used



Not made with
Latex, DEHP, or PVC

PRODUCT RANGE OVERVIEW

| Reference code | Description | Needle outer diameter | Needle length |
|----------------|---------------------------------------|-----------------------|---------------|
| HCN-R-1919-ET | Huber coreless Ndl rect 90 19 G 19 mm | 19 G | 19 mm |
| HCN-R-1925-ET | Huber coreless Ndl rect 90 19 G 25 mm | 19 G | 25 mm |
| HCN-R-2019-ET | Huber coreless Ndl rect 90 20 G 19 mm | 20 G | 19 mm |
| HCN-R-2025-ET | Huber coreless Ndl rect 90 20 G 25 mm | 20 G | 25 mm |
| HCN-R-2038-ET | Huber coreless Ndl rect 90 20 G 38 mm | 20 G | 38 mm |
| HCN-R-2219-ET | Huber coreless Ndl rect 90 22 G 19 mm | 22 G | 19 mm |
| HCN-R-2225-ET | Huber coreless Ndl rect 90 22 G 25 mm | 22 G | 25 mm |
| HCN-R-2232-ET | Huber coreless Ndl rect 90 22 G 19 mm | 22 G | 19 mm |
| HCN-R-2238-ET | Huber coreless Ndl rect 90 22G 19 mm | 22 G | 25 mm |

Needle cap lengths:

19 mm (¾") needle length: 23,6 mm needle cap

25 mm (1") needle length: 30 mm needle cap

32 mm (1 ¼") and 38 mm (1 ½") needle length: 42,7 mm needle cap

Packaging Details

| Individual unit | | | | Inner box | | | | Outer box | | | | |
|-----------------|--------|------------|----------------|-----------|--------|--------|------------|----------------|--------|--------|--------|-------------|
| L (mm) | H (mm) | Weight (g) | Quantity (pcs) | L (mm) | W (mm) | H (mm) | Weight (g) | Quantity (pcs) | L (mm) | W (mm) | H (mm) | Weight (kg) |
| 85 | 65 | 0,7 | 50 | 195 | 90 | 70 | 120 | 1000 | 414 | 387 | 211 | 3,2 |

Languages

26 languages of the IFU available:

(BG) Bulgarian, (CS) Czech, (DA) Danish, (DE) German, (EL) Greek, (EN) English, (ES) Spanish, (ET) Estonian, (FI) Finnish, (FR) French, (HR) Croatian, (HU) Hungarian, (IT) Italian, (LT) Lithuanian, (LV) Latvian, (NL) Dutch, (NO) Norwegian, (PL) Polish, (PT) Portuguese, (RO) Romanian, (RU) Russian, (SK) Slovak, (SL) Slovenian, (SR) Serbian, (SV) Swedish, (TR) Turkish.







Transport and storage conditions

Keep dry and away from sunlight.



Label details

| | Individual unit | Inner box | Outer box |
|---|-----------------|-----------|-----------|
| Product name | ✓ | ✓ | ✓ |
| Product size | ✓ | ✓ | ✓ |
| Quantity (pcs) | — | ✓ | ✓ |
| REF | ✓ | ✓ | ✓ |
| LOT | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
| STERILE EO | ✓ | ✓ | ✓ |
| CE | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |
|  | ✓ | ✓ | ✓ |

| | Individual unit | Inner box | Outer box |
|---|-----------------|-----------|-----------|
|  | ✓ | ✓ | ✓ |
| EC REP | ✓ | ✓ | ✓ |
|  | — | — | ✓ |
|  | — | — | ✓ |
|  | — | — | ✓ |
|  | — | — | ✓ |
|  | ✓ | ✓ | ✓ |
| Not made with LATEX | ✓ | ✓ | ✓ |
| Not made with DEHP | ✓ | ✓ | ✓ |

