

NIPROSOL

ACID CONCENTRATE CANISTERS



ACID CONCENTRATE CANISTERS

The ready-to-use dialysis fluid is prepared by the hemodialysis machine by mixing the acid concentrate with the bicarbonate concentrate and the purified water. The standard mixing ratios of the acid concentrates is 1+44.



NIPROSOL PORTFOLIO

	Na ⁺	Ca ⁺⁺	Mg ⁺⁺	K ⁺	Cl ⁻	CH ₃ COO ⁻	HCO ₃ ⁻	Glucose	Osmolarity
	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	mmol/l	g/l	mOsm/l
306-A (A+B)	139	1,25	0,48	1,50	103,14	2,85	37,96	1,50	294,50
307-A (A+B)	139	1,75	0,48	1,50	104,14	2,85	37,96	1,50	296,00
308-A (A+B)	139	1,50	0,48	1,50	103,64	2,85	37,96	1,00	292,48
310-A (A+B)	139	1,50	0,50	1,00	103,20	2,85	37,96	1,00	291,56
313-A (A+B)	140	1,50	0,50	1,50	108,50	3,00	34,00	1,00	294,55
314-A (A+B)	140	1,25	0,50	1,50	108,00	3,00	34,00	1,00	293,80
316-A (A+B)	140	0,63	0,25	1,50	102,44	4,00	36,81	1,00	291,18
318-A (A+B)	140	1,75	0,50	2,00	109,50	3,00	34,00	1,50	299,08
319-A (A+B)	139	1,75	0,48	3,00	105,64	2,85	37,96	1,50	299,00
320-A (A+B)	139	1,25	0,50	3,00	104,70	2,85	37,96	1,50	297,59
321-A (A+B)	140	1,75	0,50	2,00	109,50	3,00	34,00	1,00	296,30
323-A (A+B)	140	1,50	0,50	2,00	109,00	3,00	34,00	1,00	295,55
324-A (A+B)	140	1,25	0,50	2,00	108,50	3,00	34,00	1,00	294,80
326-A (A+B)	139	1,50	0,48	2,50	104,64	2,85	37,96	1,50	297,25
327-A (A+B)	140	1,50	0,50	3,00	110,00	3,00	34,00	1,00	297,55
330-A (A+B)	140	1,75	0,50	3,00	110,50	3,00	34,00	1,00	298,30
331-A (A+B)	140	1,25	0,50	3,00	109,50	3,00	34,00	1,00	296,80
332-A (A+B)	140	1,50	0,50	4,00	111,00	3,00	34,00	1,00	299,55

PACKAGING DETAILS

	Volume (L)	Article number	Canisters per box	Canisters per pallet
Stackable canisters	3,9	SOL45- [*] A4L	N/A	132
	5	SOL45- [*] A5L	N/A	132
	10,5	SOL45- [*] A11L	N/A	60
Canisters in a box	3,9	SOL45- [*] A4LBX	2	128
	5	SOL45- [*] A5LBX	2	128

^{*} acid formula

Not all products are registered and approved for sale in all countries or regions. Indications of use may also vary by country and region. Please contact your country representative for product availability and information.

SHELF-LIFE: 3 years

Nipro Renal Care is part of Nipro Corporation Japan, a leading global healthcare company established in 1954. With over 35.000 employees worldwide, Nipro serves the Medical Device, Pharmaceutical, and Pharmaceutical Packaging industries.

Nipro Renal Care is a global market leader with over 6 decades providing renal solutions for dialysis and dialysis-related treatment. We specialize in developing dialysis machines, water treatment systems, and a comprehensive portfolio of disposable medical equipment.

In order to address the needs of patients, healthcare professionals, and procurement managers alike, Nipro Renal Care is driven by innovation and patient safety to offer the highest quality products that optimize time, effort, and costs.

BECAUSE EVERY LIFE DESERVES AFFORDABLE CARE



www.nipro-group.com/en-en/our-company/our-locations

Please contact your local representative for more information.

